

## V3 SECTIONAL MATRIX SYSTEM

KEY BENEFITS & STARTER PACK CONFIGURATION



### PERFECT CONTACTS - IDEAL ANATOMY

- + Super-elastic nickel-titanium spring for optimal separation force and tight contacts
- + Wide indications for use - glass-fiber tines will not collapse into wide cavities, even when a cusp is missing
- + Tines anatomically shaped for superior matrix retention and adaptation - and minimal finishing
- + Tines' V-shape accommodates the wedge
- + V3 Tab-Matrix with pronounced marginal ridge and vertical profile for correct anatomical form
- + Wave-Wedge stackable, self-locating and self-seating for optimal seal on matrix
- + Spring tilt allows easy stacking of multiple rings
- + Ring very stable on undercuts, and in forceps during placement

#### V3 STARTER PACK INCLUDES:

- + 2 x V3 Ring Universal (green)
- + 1 x V3 Ring Narrow (yellow)
- + 25 x Wave-Wedges Small (white)
- + 25 x Wave-Wedges Medium (pink)
- + 25 x Wave-Wedges Large (purple)
- + 25 x V3 Tab-Matrices 4.5mm
- + 25 x V3 Tab-Matrices 5.5mm
- + 25 x V3 Tab-Matrices 6.5mm
- + 5 x WedgeGuard Small (white)
- + 5 x WedgeGuard Medium (pink)
- + 5 x WedgeGuard Large (purple)
- + 1 x Pin-Tweezers
- + 1 x Forceps

#### PLUS INCLUDES SUPERCURVE NON-STICK MATRIX SAMPLES:

- + 10 x SuperCurve Non-Stick Matrices 4.5mm
- + 25 x SuperCurve Non-Stick Matrices 5.5mm
- + 25 x SuperCurve Non-Stick Matrices 6.5mm

PRODUCT CODE: 403318



www.trident.com

supercurve  
NON-STICK matrix

trident  
Innovative. Simple. Smart.

Visit [www.suredental.com](http://www.suredental.com) for technique video



Call for pricing


- A) New shape with pronounced marginal ridge.
- B) An S- shaped vertical contour provides close replication of the true anatomy of the proximal surface and allows the matrix to slide past the gingival margin.
- C) V3 matrices are 30<sup>th</sup> thin
- D) Greater curvature on the horizontal plane means more of the curvature remains as the matrix wraps around the tooth.
- E) The Tab acts as a handle for easy placement.
- F) Holes in the tabs of the V3 matrix allow easy removal with pin tweezers.

The V3 Tab-Matrix takes the Tab-Matrix a step further, introducing a pronounced ridge on the occlusal margin. This curved edge produces completely natural-looking Class II composite restorations that require minimal finishing, saving time, money and effort.

Available in pkgs 50 or 100 4.5mm,5.5mm,6.5mm

## SuperCurve

"The new SuperCurve matrix from Triodent is the best matrix I have ever used. Its shape makes it incredibly easy to place, and creates perfect contours with minimal finishing. The non-stick surface makes removal a breeze. Finally, a perfect matrix!"



Small

Medium

Large

Call for pricing

## Original V-Matrix

The Tab-Matrix is anatomically accurate and easy to handle, thanks to the tab and holes which fit Triodent's Pin-Tweezers, making placement and removal simple



4.5

5.5

6.5



\$75.00

100 Pk

## Wave-Wedge

Used with the V3 Ring, the Wave-Wedge has unsurpassed sealing capability at the gingival margin. Its smart features, including a central concavity - the 'wave' - ensure it stands head and shoulders above any other wedge available.

Available in pkgs 100



1-888-277-2335